TB-CHECK-1

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
X. BIBLIOGRAFIA

1. **Grange J.M.** (1984), "The humoral immune response in Tuberculosis: its nature, biological role and diagnostic usefulness." In: Fow W., Grosset J., Styblo K. (eds) – Advances in tuberculosis research vol. 21 Karger – Basel 1:78.
2. **Krambovitis E.** (1987), "Serodiagnosis of Tuberculosis in perspective" – Serodiagn. Immunoth 1: 7-19.
3. **Daniel T.M. and Debanne S.M.** (1987), "The serodiagnosis of tuberculosis and other mycobacterial diseases by ELISA" – Am. Rev. Respir. Dis. 135:1137-1151.
4. **Toman K.** (1981) "Sensitivity, specificity and predictive value of diagnosis tests" – Bull. Int. Union Tuberc. 56:18/28.

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TB-CHECK-1

Immunochromatographic rapid test for the detection of Tuberculosis antibody in whole blood, serum or plasma samples

- Ref. 22071 -

I- PRINCIPLE

Tuberculosis remains an important socio-economical and medical problem throughout the world. In 1995, more people died of tuberculosis in the world than of any other infectious disease. According to the Center of Disease Control, the incidence of TB is expected to increase from 7.5 millions cases per year in 1995 to 11.9 millions in 2005. The case fatality rate is estimated at 55% for untreated people and 15% for the treated patients.

The commonly used diagnostic tests for the tuberculosis such as sputum examination of acid fast bacilli, culture of sputum or other fluids, the tuberculin skin test and radiological investigations, do not achieve the required diagnostic sensitivity.

Actually, the best assets for clinical and epidemiological work in the diagnosis of tuberculosis are serological tests such as haemagglutination, complement fixation, immunofluorescence, ELISA and radioimmunoassay techniques (1, 2, 3). But the highly qualified technicians and expensive material needed to perform tuberculosis tests using the above methods are a serious limitation in developing countries.

TB-CHECK-1 is a rapid qualitative screening test for the detection of human antibodies to Koch's bacillus in whole blood, serum or plasma samples. The test is able to specifically detect antibodies from patient in active phase of disease. Samples from vaccinated people will not react with TB-CHECK-1.

The method employs a unique combination of anti-human immunoglobulins dye conjugate and highly purified BCG proteins to specifically detect anti Koch antibodies. The test is mainly detecting IgG and IgA but IgM are also reacting in case of high concentration. The sensitivity has been adjusted to 350 units in accordance with a previous study (4) establishing this cut-off value as the best regarding specificity, sensitivity and predictive value for a serological tuberculosis test. As the samples flows through the absorbent device, the anti human immunoglobulins dye conjugate binds to the human IgG antibodies forming an antibody-antigen complex. This complex binds to the BCG proteins in the positive reaction zone (T) and produces a pink-rose colour band. In the absence of anti BCG antibodies, there is no line in the reaction zone (T). The reaction mixture continues flowing through the absorbent device past the reactive zone and control zone. Unbound conjugate binds to the reagents in the control zone (C), producing a pink-rose colour band, demonstrating that the reagents are functioning correctly.

II- TB-CHECK-1 KIT COMPONENTS

Each kit contains everything needed to perform 10 or 20 tests.

| | | |
|--|-------|-----|
| - TB-CHECK-1 reaction devices | 10 | 20 |
| - Diluent in a dropper bottle containing saline buffer, detergent and sodium azide (NaN ₃ < 0.1%) | 2.5mL | 5mL |
| - Disposable plastic pipettes | 10 | 20 |
| - Instruction leaflet | 1 | 1 |

III- STORAGE AND STABILITY

1- All TB-CHECK-1 kit components should be stored at room temperature (+4°C to +30°C) in the sealed pouch.

2- **Do not freeze the test kit**

3- The TB-CHECK-1 is stable until the expiry date stated on the package label.

IV- PRECAUTIONS

1- This test is designed for *in vitro* diagnostic use and professional use only.

2- Read carefully instruction leaflet before using this test.

3- Do not use beyond the expiry date which appears on the package label.

4- Do not use a test from a damaged protective wrapper.

5- Handle all specimens as if they contained infectious agents. When the assay procedure is completed, dispose of specimens carefully after autoclaving them for at least one hour. Alternatively, they can be treated with 0.5 to 1% solution of sodium hypochlorite for one hour before disposal.

6- Wear protective clothing such as laboratory coats and disposable gloves while assaying samples.

7- Do not eat, drink or smoke in the area where specimens and kit reagents are handled.

8- Avoid any contact between hands and eyes or nose during specimen collection and testing.

V- SPECIMEN COLLECTION AND PREPARATION

Serum, plasma (lithium or ammonium heparinate, EDTA) or whole blood

The specimen should be collected under the standard laboratory conditions (aseptically in such a way as to avoid hemolysis).

When the test is performed with whole blood, fresh samples should be used (< 4hours).

If the test is to be run within 48 hours after collection, the specimen should be stored in the refrigerator (+2°C to +8°C). If testing is delayed more than 48 hours, the specimen should be frozen. The frozen specimen must be completely thawed,

thoroughly mixed and brought to room temperature prior to testing. Avoid repeated freezing and thawing.

In case of cloudiness, high viscosity or presence of particulate matter into the serum specimen, it should be diluted with equal volume (V/V) of diluting buffer (not provided but available upon request) before testing.

VI- ASSAY PROCEDURE


- 1- Allow samples and TB-CHECK-1 test devices to come to room temperature prior to testing.
- 2- Remove the reaction device from its protective wrapper by tearing along the split.
- 3- Label device with the patient's name or control number.
- 4- Fill the serum dropper with specimens (serum or plasma) and by holding it vertically, dispense one drop (25 µL) into sample well (▷). If whole blood is used, dispense two drops (50 µL) into the sample well (▷) **and wait for the whole blood sample to be completely absorbed before adding diluent.**
- 5- Add exactly 5 full drops of diluent (200 µL) in the sample well (▷).
- 6- Read the results after 10 to 15 minutes.

DO NOT INTERPRETE THE TEST AFTER 15 MINUTES.

VII- READING TEST RESULTS


A. Negative

Only one coloured band shows in the control zone (C).

B. Positive

In addition to the control band, a clearly distinguishable band also shows in the test zone (T) indicating a positive result.

C. Inconclusive

If there is no distinct coloured band visible in the control zone (C), the test is inconclusive. In this case, it is recommended to repeat the test.

VIII- PERFORMANCES CHARACTERISTICS

a) Sensitivity

A panel of 23 positive and negative sera, preassayed with an ELISA method for IgG, IgA and IgM, has been tested with TB-CHECK-1. The results are summarized in Table I and show a sensitivity of around 250 U/mL both for IgG and IgA. The final sensitivity test has been further adjusted to 350 U/mL in order to optimize the sensitivity, the specificity and the predictive value of the test (4).

| Serum | ELISA | | | TB-CHECK-1 (15 min) |
|---------|------------|------|------------|---------------------|
| | IgG (U/mL) | IgA | IgM (U/mL) | |
| 830 | neg. | neg. | neg. | neg. |
| 866 | neg. | neg. | neg. | neg. |
| 868 | neg. | neg. | neg. | neg. |
| 869 | neg. | neg. | neg. | neg. |
| 870 | neg. | neg. | neg. | neg. |
| 871 | neg. | neg. | neg. | neg. |
| 876 | neg. | neg. | neg. | neg. |
| 881 | neg. | neg. | neg. | neg. |
| 883 | neg. | neg. | neg. | neg. |
| 910 | neg. | neg. | neg. | neg. |
| 919 | neg. | neg. | neg. | neg. |
| 921 | neg. | neg. | neg. | neg. |
| HM | neg. | neg. | neg. | neg. |
| 832 | 130 | neg. | weak+ | pos. |
| 873 | 258 | neg. | neg. | pos. |
| 834 | 430 | neg. | weak+ | pos. |
| 879 | 525 | 200 | neg. | pos. |
| 849 | neg. | 290 | pos. | pos. |
| 878 | weak+ | 435 | pos. | pos. |
| 828 | neg. | 240 | pos. | pos. |
| 849 | neg. | 290 | pos. | pos. |
| M 10/11 | 700 | 350 | weak+ | pos. |
| M 24/10 | 1000 | 1000 | pos. | pos. |
| 874 | 500 | 1600 | pos. | pos. |

Table I : Comparison of TB-CHECK-1 and Elisa method

Another study was performed on 391 clinically characterized people. Results are summarized in Table II.

| Clinical status | TB-CHECK-1 results | | | |
|--|--------------------|----------|----------|------------|
| | n. | Positive | Negative | Positive % |
| 1. Active pulmonary tuberculosis | 130 | 113 | 17 | 87 |
| 2. Active non tuberculosis mycobacteriosis (NTM) | 24 | 20 | 4 | 83 |
| 3. Old tuberculosis | 15 | 4 | 11 | 27 |
| 4. Other respiratory diseases | 79 | 4 | 75 | 5 |
| 5. Other diseases | 60 | 2 | 58 | 3 |
| 6. Healthy subjects | 83 | 7 | 76 | 8 |
| Total | 391 | | | |

Table II : TB-CHECK-1 evaluation results

The sensitivity of the test is shown to be 87% (113/130) for pulmonary tuberculosis, 83% (20/24) non tuberculosis mycobacteriosis and 86% (133/154) overall.

The Tables III, IV, V and VI show the performances of TB-CHECK-1 in comparison with different detection methods (smear test, sputum culture or PPD) and according to clinical recordings. Only people belonging to the groups (1) or (2) (active mycobacteriosis) were selected for these studies.

| Smear test results | TB-CHECK-1 Results | | | |
|--------------------|--------------------|------------|-----------|-----------------|
| | n. | Positive | Negative | Sensitivity (%) |
| G0 (negative) | 23 | 18 | 5 | 78 |
| G1-G3 | 66 | 54 | 12 | 82 |
| G4-G6 | 46 | 43 | 3 | 93 |
| G7-G10 | 19 | 18 | 1 | 95 |
| Total | 154 | 133 | 21 | 86.4% |

Table III : Comparison of the results of

| Sputum culture results | TB-CHECK-1 Results | | | |
|------------------------|--------------------|------------|-----------|-----------------|
| | n. | Positive | Negative | Sensitivity (%) |
| - | 6 | 5 | 1 | 83 |
| + | 75 | 59 | 16 | 79 |
| 2+ | 35 | 32 | 3 | 91 |
| 3+ | 23 | 22 | 1 | 96 |
| 4+ | 11 | 11 | 0 | 100 |
| Total | 150 | 129 | 21 | 86% |

Table IV : Comparison of the results of TB-CHECK-1 and sputum culture methods

| PPD reaction results | TB-CHECK-1 Results | | | |
|----------------------|--------------------|-----------|-----------|-----------------|
| | n. | Positive | Negative | Sensitivity (%) |
| 0.5-5 mm | 5 | 4 | 1 | 80 |
| 6-10 mm | 7 | 6 | 1 | 86 |
| 11-20 mm | 33 | 24 | 9 | 73 |
| 21-40 mm | 29 | 25 | 4 | 86 |
| 41- mm | 10 | 9 | 1 | 90 |
| Total | 84 | 68 | 16 | 81% |

Table V : Comparison of the results of TB-CHECK-1 and PPD reaction methods

| Clinical recordings | | TB-CHECK-1 RESULTS | | | |
|----------------------|----------------------|--------------------|------------|-----------|-----------------|
| | | n. | Positive | Negative | Sensitivity (%) |
| First attack : | - Less than 3 months | 64 | 53 | 11 | 83 |
| | - More than 3 months | 29 | 27 | 2 | 93 |
| Relapse : | - Less than 3 months | 15 | 13 | 2 | 87 |
| | - More than 3 months | 10 | 9 | 1 | 90 |
| Chronic tuberculosis | | 12 | 11 | 1 | 92 |
| Total | | 130 | 113 | 17 | 87% |

Table VI : TB-CHECK-1 results according to clinical recordings

b) Specificity

From results in Table II, the specificity of the test was calculated as follows :

$$\begin{aligned}
 \text{Specificity} &= \frac{\text{negatives with TB-CHECK-1}}{\text{negatives by clinical status}} \times 100 \\
 &= \frac{(11+75+58+76)}{(15+79+60+83)} \times 100 \\
 &= \frac{220}{237} \times 100 = 92.8\%
 \end{aligned}$$

c) Overall agreement

From results summarized in Table II, the overall agreement of the test regarding the clinical status was calculated as follows :

$$\begin{aligned}
 \text{Correlation} &= \frac{\text{true positives} + \text{true negatives}}{\text{total cases}} \times 100 \\
 &= \frac{(113+20) + (11+75+58+76)}{391} \times 100 \\
 &= \frac{353}{391} \times 100 = 90\%
 \end{aligned}$$

d) Cross-reaction

As shown in table II, TB-CHECK-1 is highly cross-reacting with other mycobacteria infections (groupe 2 : active non tuberculosis mycobacteriosis).

e) Reproducibility

Intra-assay

Within run precision was determined by using 3 replicates of four positive (different antibodies level) specimens and one negative specimen. The negative and positives values were correctly identified 100% of time.

Inter-assay


Between run precision was determined by using the same four positive and negative specimens as for intra –assay with five different lots of reaction device. The negative and positive values were correctly identified 100% of time.

IX- LIMITATIONS OF THE TEST

- 1- As for any *in vitro* diagnostic procedure, the physician should evaluate data obtained by the use of the test in light of other clinical information.
- 2- Very early TB infections could lead to false negative results, due to the low concentration of anti-TB antibodies in the serum, plasma or whole blood samples.
- 3- TB-CHECK-1 is showing high cross-reaction with other mycobacteria and such cases should be identified before testing.
- 4- A positive result does not exclude the presence of other pathogens.

X- BIBLIOGRAPHY

- 1- **Grange J.M.** (1984), "The humoral immune response in Tuberculosis : its nature, biological role and diagnostic usefulness." In : Fow W., Grosset J., Styblo K. (eds) - Advances in tuberculosis research vol.21 Karger - Basel 1:78.
- 2- **Krambovitis E.** (1987), "Serodiagnosis of Tuberculosis in perspective" - Serodiagn. Immunoth 1 : 7-19.
- 3- **Daniel T.M. and Debanne S.M.** (1987), "The serodiagnosis of tuberculosis and other mycobacterial diseases by ELISA" - Am. Rev. Respir. Dis. 135 : 1137-1151.
- 4- **Toman K.** (1981) "Sensitivity, specificity and predictive value of diagnosis tests" - Bull. Int. Union Tuberc. 56 : 18/28.

Temperature limitations

Consult operating instructions

Do not re-use

For *in vitro* diagnostic use

Manufactured by VEDALAB - France