
Test rapido Anticorpo *H.pylori* (Siero/Plasma)

Foglietto Illustrativo

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【BIBLIOGRAFIA】

- Marshall, BJ, McGeachie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. (1985), 149: 439-44.
- Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. (1990), 322: 909-16.
- Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. (1987), 82(4): 292-96.
- Loffeld, RJLF, et al. Usefulness of several commercial enzyme-linked immunoassays for detection of Helicobacter pylori infection in clinical medicine. Euro. J. Gastroen. Hepa. (1993) 5:333-37.
- Cutler, AF, et al. Accuracy of invasive and non-invasive tests to diagnose Helicobacter pylori infection. Gastroenterology.(1995), 109: 136-141.
- Ansorg, R, Von Recklinghausen, G, Pomarius, R and Schmid, EN. Evaluation of techniques for isolation, subcultivation and preservation of Helicobacter pylori. J. Clin. Micro. (1991), 29:51 53.
- Pronovost, AP, Rose, SL, Pawlak, J, Robin, H and Schneider, R. Evaluation of a new immunodiagnostic assay for Helicobacter pylori antibody detection: Correlation with histopathological and microbiological results. J. Clin. Micro. (1994), 32: 46-50.
- Megraud, F, Bassens-Rabbe, MP, Denis, F, Belbourni, A and Hoa, DQ. Seroprevalence of Campylobacter pylori infection in various populations. J. Clin. Micro. (1989), 27: 1870-3.

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ACRO Biotech, Inc.

9500 Seventh Street,
Unit M, Rancho Cucamonga,
CA 91730, U.S.A.

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For professional in vitro diagnostic use only.

【INTENDED USE】

The H.pylori Antibody Rapid Test Cassette (Serum/Plasma) is a rapid chromatographic immunoassay for the qualitative detection of Antibody to H. pylori in serum or plasma.

【SUMMARY】

H.pylori is a small, spiral-shaped bacterium that lives in the surface of the stomach and duodenum. It is implicated in the etiology of a variety of gastrointestinal diseases, including duodenal and gastric ulcer, non-ulcer dyspepsia and active and chronic gastritis.^{1,2} Both invasive and non-invasive methods are used to diagnose H.pylori infection in patients with symptoms of gastrointestinal disease. Specimen-dependent and costly invasive diagnostic methods include gastric or duodenal biopsy followed by urease testing (presumptive), culture, and/or histologic staining.³ Non-invasive techniques include the urea breath test, which requires expensive laboratory equipment and moderate radiation exposure, and serological methods.^{4,5} Individuals infected with H.pylori develop antibodies which correlate strongly with histological confirmed H.pylori infection.^{6,7,8}

The H.pylori Antibody Rapid Test Cassette (Serum/Plasma) is a simple test that utilizes a combination of H.pylori antigen coated particles and anti-human IgG to qualitatively and selectively detect H.pylori antibodies in serum or plasma.

【PRINCIPLE】

The H.pylori Antibody Rapid Test Cassette (Serum/Plasma) is a qualitative membrane based immunoassay for the detection of H.pylori antibodies in serum or plasma. In this test procedure, anti-human IgG is immobilized in the test line region of the test. After specimen is added to the specimen well of the cassette, it reacts with H.pylori antigen coated particles in the test. This mixture migrates chromatographically along the length of the test and interacts with the immobilized anti-human IgG. If the specimen contains H.pylori antibodies, a colored line will appear in the test line region indicating a positive result. If the specimen does not contain H. pylori antibodies, a colored line will not appear in this region indicating a negative result. To serve as a procedural control, a colored line will always appear in the control line region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

【REAGENTS】

The test cassette contains H. pylori antigen coated particles and anti-human IgG coated on the membrane.

【PRECAUTIONS】

Please read all the information in this package insert before performing the test.

- For professional in vitro diagnostic use only. Do not use after the expiration date.
- Do not eat, drink or smoke in the area where the specimen or kits are handled.
- Handle all the specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
- Humidity and temperature can adversely affect results.

【STORAGE AND STABILITY】

Store as packaged at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

【SPECIMEN COLLECTION AND PREPARATION】

- The H.pylori Antibody Rapid Test Cassette (Serum/Plasma) can be performed using serum or plasma.
- Separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear, non-hemolyzed specimens.
- Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8°C for up to 3 days. For long term storage, specimens should be kept below -20°C.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- If specimens are to be shipped, they should be packed in compliance with federal regulations covering the transportation of etiologic agents.

【MATERIALS】

Materials provided

- Test cassettes
- Droppers
- Package insert

Materials required but not provided


- Specimen collection containers
- Centrifuge
- Timer

【DIRECTIONS FOR USE】

Allow test cassette, serum or plasma specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.

- Remove the test cassette from the sealed foil pouch and use it as soon as possible. Best results will be obtained if the assay is performed immediately after opening the foil pouch.

- Hold the dropper vertically and transfer **3 drops of serum or plasma (approximately 75 µL)** to the specimen well of test cassette and start the timer. Avoid trapping air bubbles in the specimen well. See illustration below.
- Wait for the colored line is appeared. The result should be read **at 10minutes**. Do not interpret the result after **20 minutes**.

【INTERPRETATION OF RESULTS】

(Please refer to the illustration above)

POSITIVE: * **Two distinct colored lines appear.** One colored line should be in the control region (C) and another colored line should be in the test region (T).

***NOTE:** The intensity of the color in the test line region (T) will vary depending on the concentration of H.pylori antibody present in the specimen. Therefore, any shade of color in the test region (T) should be considered positive.

NEGATIVE: **One colored line appears in the control region (C).** No apparent colored line appears in the test region (T).

INVALID: **Control line fails to appear.** Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

【QUALITY CONTROL】

Internal procedural controls are included in the test. A red line appearing in the control region (C) is an internal positive procedural control. It confirms sufficient specimen volume and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

【LIMITATIONS】

- The H.pylori Antibody Rapid Test Cassette (Serum/Plasma) is for in vitro diagnostic use only. This test should be used for the detection of H.pylori antibody in serum or plasma specimen. Neither the quantitative value nor the rate of increase in H.pylori antibody concentration can be determined by this qualitative test.
- The H.pylori Antibody Rapid Test Cassette (Serum/Plasma) will only indicate the presence of H.pylori antibodies in the specimen and should not be used as the sole criteria for the diagnosis of H.pylori infection.
- As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of H.pylori infection.

【EXPECTED VALUES】

The H.pylori Antibody Rapid Test Cassette (Serum/Plasma) has been compared with a leading commercial H.pylori antibody EIA test. The correlation between these two systems is over 97%.

【PERFORMANCE CHARACTERISTICS】

Sensitivity and Specificity

The H.pylori Antibody Rapid Test Cassette (Serum/Plasma) has been evaluated with serum and plasma specimens obtained from a population of symptomatic and asymptomatic individuals who presented for endoscopic examination.

Method	ELISA			Total Results
	Results	Positive	Negative	
H.pylori Antibody Rapid Test Cassette(Serum/Plasma)	Positive	171	8	179
	Negative	0	102	102
Total Results		171	110	281

Relative Sensitivity: >99.9% (95%CI*: 98.3%-100%) *Confidence Interval

Relative Specificity: 92.7% (95%CI*: 86.2%-96.8%)

Accuracy: 97.2% (95%CI*: 94.5%-98.8%)

Precision

Intra-Assay

Within-run precision has been determined by using 10 replicates of 4 specimens: a negative, a low positive, a medium positive and a high positive. The negative, low positive, medium positive and high positive values were correctly identified >99% of the time.

Inter-Assay

Between-run precision has been determined by 10 independent assays on the same 4 specimens: a negative, a low positive, a medium positive and a high positive. Three different lots of the H. pylori Test Cassette (Serum/Plasma) have been tested using negative, low positive medium positive and high positive specimens. The specimens were correctly identified >99% of the time.

Cross-reactivity

Sera containing known amounts of antibodies to H.pylori have been tested with Hepatitis A, B, C, E, HIV and Syphilis. No cross-reactivity was observed, indicating that the H.pylori Test Cassette (Serum/Plasma) has a high degree of specificity for human antibodies to H.pylori.

Interfering Substances

The H.pylori Antibody Rapid Test Cassette (Serum/Plasma) has been tested for possible interference from visibly hemolyzed and lipemic specimens, as well as specimens containing high bilirubin levels. In addition, no interference was observed in specimens containing up to 1,000 mg/dL hemoglobin, up to 1,000 mg/dL bilirubin, and up to 2,000 mg/dL human serum albumin.

【BIBLIOGRAPHY】

- Marshall, BJ, McGeachie, DB, Rogers, PAR and Glancy, RG. Pyloric Campylobacter infection and gastroduodenal disease. Med. J. Australia. 149: 439-44; 1985.
- Soll, AH. Pathogenesis of peptic ulcer and implications for therapy. New England J. Med. 322:909-16; 1990.
- Hazell, SL, et al. Campylobacter pyloridis and gastritis I: Detection of urease as a marker of bacterial colonization and gastritis. Amer. J. Gastroenterology. 82(4): 292-96; 1987.
- Loffeld, RJLF, et al. Usefulness of several commercial enzyme-linked immunoassays for detection of Helicobacter pylori infection in clinical medicine. Euro. J. Gastroen. Hepa. 5:333-37; 1993.
- Cutler, AF, et al. Accuracy of invasive and non-invasive tests to diagnose Helicobacter pylori infection. Gastroenterology. 109: 136-141; 1995.
- Ansorg, R, Von Recklinghausen, G, Pomarius, R and Schmid, EN. Evaluation of techniques for isolation, subcultivation and preservation of Helicobacter pylori. J. Clin. Micro. 29:51-53; 1991.
- Pronovost, AP, Rose, SL, Pawlak, J, Robin, H and Schneider, R. Evaluation of a new immunodiagnostic assay for Helicobacter pylori antibody detection: Correlation with histopathological and microbiological results. J. Clin. Micro. (1994), 32: 46-50.
- Megraud, F, Bassens-Rabbe, MP, Denis, F, Belbourni, A and Hoa, DQ. Seroprevalence of Campylobacter pylori infection in various populations. J. Clin. Micro. 27: 1870-3; 1989.

Index of Symbols

	Attention, see instructions for use		Tests per kit		Authorized Representative
	For in vitro diagnostic use only		Use by		Do not reuse
	Store between 2-30°C		Lot Number		Catalog #
	Do not use if package is damaged				

ACRO Biotech, Inc.
9500 Seventh Street,
Unit M, Rancho Cucamonga,
CA 91730, U.S.A.

MedNet GmbH
Borkstrasse 10
48163 Muenster
Germany

Number: 145438400
Effective date: 2016-12-05