PSA-CHECK-1

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X. BIBLIOGRAFIA

1–**Bagshawe, K.D.** 1993. Tumor markers. Br J. Cancer 48 : 167-175.

2–**Berg, W., Linder, Ch., Eschholz, G., Link, St., Schubert, J.** 1999. Possibility of Improving the Acceptance Rate of Early Detection Testing for Prostate Cancer with a OneStep Test for ProstateSpecific Antigen in Whole Blood. Urol. Int. 63 : 102–106.

3– **Irani, J. Millet, C., Levillain, P., Doré, B., Begon, F., Aubert, J.** 1997. Serurtourinary prostate specific antigen ratio : its impact in distinguishing prostate cancer when serum prostate specific antigen level is 4 to 10 ng/ml. J. of Urology. 157 : 185-188.

4–**Kuriyama, M, MC Wang, CL Lee, LD Papsidero, C.S. Killian, H. Inaji, N.H. Slack, T. Nishiura, GP. Murphy and T.M. Chu.** 1981. Use of human prostate specific antigen in monitoring cancer. Cancer res. 41: 3874-3876.


5–**Liedtke R.L. and JD Batjer.** 1984. Measurement of prostate specific antigen by radioimmunoassay. Clin. Chem. 30 : 649-652.

6–**Papsidero, LD, G.A. Crighan, M.C. Wang, M. Kuriyama, E.A. Johnson, L.A. Valenzuela and T.M. Chu.** 1983. Monoclonal antibody to human prostate antigen. Hybridoma 2 : 139-147.

7–**Wang, M.C., M. Kuriyama, L.D. Papsidero, R.M. Loor, L.A. Valenzuela, G.P. Murphy and T.M. Chu.** 1982. Prostate antigen of human cancer patients, P. 179-197.

In H. Busch and L.C. Yeoman (ed). Methods in cancer research. Academic Press Inc, New York.

8–**Piironen T., Pettersson K., Suonpää M., Stenman U.H., Oesterling J. E., Lövgren T. and Lilja H.** 1996. In vitro stability of free ProstateSpecific Antigen (PSA) and ProstateSpecific Antigen (PSA) complexed to α1antichymotrypsin in blood samples. Urology 48 (6A) : 81-87.

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PSA-CHECK-1

Qualitative rapid test for the detection of Prostate-Specific Antigen

- Whole Blood, plasma or serum –
Ref. 8051A

I- PRINCIPLE

Prostate-specific antigen (PSA) is an intracellular glycoprotein (M.W. 34000 Dalton) synthesized only by the prostate gland and in seminal plasma. PSA, a normal constituent of prostate tissue, is present in benign hyperplastic and malignant prostatic tissue, in metastatic prostatic carcinoma, and in prostatic fluid and seminal plasma. PSA is elevated in serum of prostate cancer patients due to the release of the antigen into circulation.

The PSA-CHECK-1 test is a rapid qualitative assay for the detection of Human PSA (free form or complexed with α -1-antichymotrypsin) in plasma, serum or whole blood.

The method employs an unique combination of monoclonal-dye (mouse) conjugate and monoclonal (mouse) solid phase antibodies to selectively identify PSA in the test samples with a high degree of sensitivity.

As the test sample flows through the absorbent device, the labelled antibody-dye conjugate binds to the PSA forming an antibody antigen complex. This complex binds to the Anti-PSA antibody in the positive reaction zone (B) and produces a pink-rose colour band when PSA concentration is higher than 4 ng/mL. In the absence of PSA, there is no line in the positive reaction zone. The reaction mixture continues flowing through the absorbent device past the reaction zone and control zone. Unbound conjugate binds to the reagents in the control zone producing a pink-rose colour band, demonstrating that the reagents are functioning correctly.

II- PSA-CHECK-1 KIT COMPONENTS

Each kit contains everything needed to perform 10 or 20 tests:

- PSA-CHECK-1 reaction devices	10	20
- Diluent in a dropper bottle containing saline buffer, detergent and sodium azide ($\text{NaN}_3 < 0.1\%$)	2.5 mL	5 mL
- Disposable plastic pipettes	10	20
- Instruction leaflet	1	1

III- STORAGE AND STABILITY

1- All PSA-CHECK-1 kit components should be stored at room temperature (+4°C to +30°C).

2- Do not freeze the test kit.

3- PSA-CHECK-1 is stable until the expiry date stated on the package label.

IV- PRECAUTIONS

- 1- For *in vitro* diagnostic use and professional use only.
- 2- Handle all specimens as if they contain infectious agents. When the assay procedure is completed, dispose of specimens carefully after autoclaving them for at least one hour. Alternatively, they can be treated with 0.5 to 1% solution of Sodium hypochlorite for one hour before disposal.
- 3- Wear protective clothing such as laboratory coats and disposable gloves while assaying samples. Avoid any contact between hands and eyes or nose during specimens collection and testing.
- 4- Do not eat, drink or smoke in the area where specimens and kit reagents are handled.
- 5- Read carefully the instruction notice before using this test.
- 6- Do not use beyond the expiration date which appears in the package label.
- 7- Do not use a test from a damaged protective wrapper.

V- SPECIMEN COLLECTION AND PREPARATION

Serum, plasma (lithium or ammonium heparinate, EDTA) or whole blood

The specimen should be collected under the standard laboratory conditions (aseptically in such a way as to avoid hemolysis).

Each specimen should be treated as if potentially infectious

Whole blood samples should be tested immediately.

If the test is to be run within 48 hours after collection the specimen should be stored (8) in the refrigerator (+2° to +8°C). If testing is delayed more than 48 hours, the specimen should be frozen (8). The frozen specimen must be completely thawed, thoroughly mixed and brought to room temperature prior to testing. Avoid repeated freezing and thawing.

In the case of cloudiness, high viscosity or presence of particulate matter into the serum specimen, it should be diluted with equal volume (V/V) of diluting buffer (not provided but available upon request) before testing.


VI- ASSAY PROCEDURE

- 1- Allow samples and PSA-CHECK-1 test devices to come to room temperature prior to testing.
- 2- Remove the “reaction device” from its protective wrapper by tearing along the split.
- 3- Label device with the patient’s name or control number.

4– Fill the serum dropper with specimens (serum, plasma or whole blood) and by holding it vertically, dispense exactly one drop (25 µL) into sample well (A). If whole blood is used, dispense exactly 2 drops (50 µL) into the sample well (A) **and wait for the whole blood sample to be completely absorbed before adding diluent.**

Warning !

In order to obtain better results, the sample (serum, plasma or whole blood) must be added at the right side of the sample well (A).


5– Add exactly 5 or 6 full drops of diluent (200 µL) into the sample well (A).

6– Read the results at 10 minutes. Do not interpret after 15 minutes.

VII– READING TEST RESULTS


Negative

Only one coloured band appears in the control window (C). No apparent band in the test region.

Positive

In addition to the control band, a clearly distinguishable band also appears in the test window (B).

Inconclusive :

If there is no distinct colour band visible in the test window and test control window, the test is inconclusive. It is recommended in this case to repeat the test.

VIII– PERFORMANCE CHARACTERISTICS

A) ACCURACY

A study was performed on 190 samples (whole blood – table 1 and serum – table 2) pretested with an EIA method. The results show a good correlation between whole blood and serum samples and also a good overall correlation with EIA results.

Table 1 : Comparison of Vedalab PSA rapid test versus EIA method with whole blood samples

PSA EIA		VEDALAB (%)	
		Reading time	
(ng/mL)	Number	5 min	10 min
0 – < 2	24	100	100
2 – < 3	29	96.6	82.8
3 – < 4	21	95.2	66.7
4 – < 5	24	41.7	70.8
5 – < 6	21	38.1	88.7
6 – < 8	27	61.5	96.3
8 – < 10	19	73.7	100
> 10	25	92.0	100

Table 2 : Comparison of Vedalab PSA rapid test versus EIA method with serum samples

PSA EIA		VEDALAB (%)	
		Reading time	
(ng/mL)	Number	5 min	10 min
0 – < 2	24	100	100
2 – < 3	29	100	82.8
3 – < 4	21	95.2	57.1
4 – < 5	24	4.2	79.2
5 – < 6	21	4.8	90.5
6 – < 8	27	33.3	88.9
8 – < 10	19	57.9	94.7
> 10	25	76.0	100

B) INTERFERENCES

Potentially interfering substances were added to a negative human serum previously tested using the PSA test. In each case, no interference with the expected Vedalab rapid PSA-check-1 test results were observed.

Kallikrein	1 µg/mL
Prostatic acid phosphatase	1 U/L
Rheumatoid factor	36 U/mL

C) PRECISION

Intra-assay

Within run precision was determined by using 10 replicates of three specimens containing 0, 5 and 10 ng/mL of PSA. The negative and positives values were correctly identified 100% of time.

Inter-assay

Between run precision was determined by using the same three specimens of 0, 5 and 10 ng/mL of PSA in 10 independent assays and with three different lots of reaction device over a 6 months period. Again, the negative and positive values were correctly identified 100% of time.

D) SENSITIVITY

PSA-CHECK-1, Rapid Test is capable of detecting PSA levels higher than 4 ng/L according to INTERNATIONAL PSA certified reference material (CRM 613 N°1004 from the Community Bureau of Reference – Belgium 1998).

However, false positive and false negative results can be found especially in the specimens around the cut-off value of 4.0 ng/mL (cf. table 1).

Serum from healthy men will normally show undetectable levels of PSA when tested on Vedalab's PSA-CHECK-1 Rapid Test.

E) HOOK EFFECT

Specimens containing very high levels of PSA (up to 10 µg/mL) when tested consistently gave positive results.

IX- LIMITATIONS

1- Some serum specimens with high rheumatoid factor concentration (RF) may yield non specific positive results during testing. Such cases should be discriminated before testing.

2- The test is designed to eliminate the potential interference of human antibodies to murine IgG (HAMA). However high level of HAMA could give falsely positive results.

3- A study (3) performed among 130 patients (53 to 95 years old) showed that a prostate cancer was confirmed for 56 patients having a serum PSA level within 1.5 up to 1515 ng/mL. BPH (benign prostatic hyperplasia) was diagnosed for 74 subjects having a serum PSA level within 0.1 and 34 ng/mL. It is advisable to determine your own values of reference for PSA level.

4- As it is true with any diagnostic procedure, the physician should evaluate data obtained by the use of this kit in light of other clinical information.

X- BIBLIOGRAPHY

1- **Bagshawe, K.D.** 1993. Tumor markers. Br J. Cancer 48 : 167-175.

2- **Berg, W., Linder, Ch., Eschholz, G., Link, St. , Schubert, J.** 1999. Possibility of Improving the Acceptance Rate of Early Detection Testing for Prostate Cancer with a One-Step Test for Prostate-Specific Antigen in Whole Blood. Urol. Int. 63 : 102- 106.

3- **Irani, J. Millet, C., Levillain, P., Doré, B., Begon, F., Aubert, J.** 1997. Serum-to-urinary prostate specific antigen ratio : its impact in distinguishing prostate cancer when serum prostate specific antigen level is 4 to 10 ng/ml. J. of Urology. 157 : 185-188.


4- **Kuriyama, M, MC Wang, CL Lee, LD Papsidero, C.S. Killian, H. Inaji, N.H. Slack, T. Nishiura, GP. Murphy and T.M. Chu.** 1981. Use of human prostate specific antigen in monitoring cancer. Cancer res. 41: 3874-3876.

5- **Liedtke R.L. and JD Batjer.** 1984. Measurement of prostate specific antigen by radioimmunoassay. Clin. Chem. 30 : 649-652.

6- **Papsidero, LD, G.A. Crighan, M.C. Wang, M. Kuriyama, E.A. Johnson, L.A. Valenzuela and T.M. Chu.** 1983. Monoclonal antibody to human prostate antigen. Hybridoma 2 : 139-147.

7- **Wang, M.C., M. Kuriyama, L.D. Papsidero, R.M. Looor, L.A. Valenzuela, G.P. Murphy and T.M. Chu.** 1982. Prostate antigen of human cancer patients, P. 179-197. In H. Busch and L.C. Yeoman (ed). Methods in cancer research. Academic Press Inc, New York.

8- **Piironen T., Pettersson K., Suonpää M., Stenman U.-H., Oesterling J. E., Lövgren T. and Lilja H.** 1996. In vitro stability of free Prostate-Specific Antigen (PSA) and Prostate-Specific Antigen (PSA) complexed to α_1 -antichymotrypsin in blood samples. Urology 48 (6A) : 81-87.

Temperature limitation

Consult operating instructions

Do not re-use

For *in vitro* diagnostic use

Manufactured by VEDALAB - France