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
The pipette can be used to withdraw blood from any kind of test-tube. Direct suction system which allows the lab-technician, by pulling the internal piston upwards, to let the blood column reach the "0" level, with automatic stop. The pipette is graduated from "0" to "180", therefore, thanks to the transparent back-ground, the reading is immediate. The rack can hold ten pipettes in pre-numbered position and ensures every pipette a perfectly vertical position.

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USE

1. Insert the pipette (BVS012) down to the bottom of the test tube containing blood and solution of sodium citrate
2. Aspirate the blood by pulling the piston up to the stop sign or the "0" and pulling the excess part of the piston
3. Fit the tube vertically in the holder support tubes (BVS011)
4. After an hour you can clearly read the erythrocytes falling index on the scale.

To make quick and safe operation can arrange the tubes directly into the support and then proceed to suction;

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IN VITRO DIAGNOSTIC MEDICAL DEVICE CE
for professional use only

- Use only for laboratory analysis
- Use gloves: contamination risk if contents leaks
- Keep out of flame or heat sources which might damage the product
- Do not use after expiry date
- Do not vary the intended purpose of the product
- Keep out of reach of children
- Store in dry place, temperature range: min. - 10°C max 50°C
- Disposal: use appropriate personal protective equipment and act according to applicable regulations

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DESTINATION USE

Disposable labware intended for use in the laboratory for aspiration of blood and subsequent detection of the erythrocyte sedimentation rate. **FOR IN VITRO USE ONLY.** The use and the manipulation of the devices **are reserved to the medical and technical staff** with the necessary professional qualification.

CLASSIFICATION

In order to the Directive 98/79/CE concerning medical-diagnostic in vitro devices receipt in Italy with D.L. n. 332 G.U.R.I. n. 269 dated 17th November 2000, Biosigma SpA defines all the components of the above-state family as medical-diagnostic in vitro devices and classifies it as "Other medical-diagnostic in vitro devices", in order to make conformable the attachment III of the directive.

The Declaration of Conformity in accordance to the **Directive 98/79/CE** is available, on request, to the Quality Department Biosigma SpA

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DEFINIZIONE DEL FABBRICANTE / MANUFACTURER DEFINITION**BIO SIGMA SpA**