Biosigma SpA

a Dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) |

Fax ++39 0426 302228 | SMS ++39 348 4077376 |

E-mail info@biosigmaeu.com | http://www.biosigma.com |

M521D

06/10/2020

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
FOR THE STANDARDIZATION OF MICROSCOPIC URINALYSIS


FAST - READ 102® system improves the standardization of microscopic urinalysis and provides precision and reproducibility such as to guarantee constant readings which will not be influenced by variations of techniques among different operators. FAST - READ 102® is made of a slide, protected by an optically transparent film, with 10 independent chambers containing a standard volume of 7 µL. After dispensing the sample on the slide application area by means of a capillary mechanism, the sediment is homogeneously distributed in the reading chamber. Each chamber is fitted with its own system for the collection of excess urine to prevent any possible contamination. Furthermore FAST - READ 102® S slide allows easy identification of the sample by using the numbers printed on the sides.

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
After completing the preparation of the urinary sediment by normal centrifugation,

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Use the pipette to mix the urinary sediment by repeatedly pressing the bulb; then

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*HOW TO USE THE NUMBERING SYSTEM
Example: For the identification of samples*

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

insert a pipette into the test tube and dispense one drop of STAIN.

dispense a drop of sample onto the appropriate area of the slide.

from number 131 to 140 mark number 1 on the left side and number 3 on the right side of the slide.

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GB - FAST READ 102[®] CELL COUNT METHOD

Fast Read 102[®] is a plastic device composed of 10 counting chambers. Each chamber contains a GRID with 10 SQUARES, subdivided into 16 small squares (called SECTORS) .

After filling the counting chamber with the sample, proceed to the counting of cells distributed in N squares.

Considering that the grid contains 10 squares and each square has a dimension of 1 x 1 mm, a depth of 0.1 mm and a volume of 0.1 ul, the formula for determining the concentration of cells (cells / ml) is:

$$[\text{Cells/ml}] = \frac{(\sum \text{ cells counted in N squares})}{N} \times \text{dilution factor} \times 10^4$$

10^4 = conversion from 0.1ul to 1 ml

Attention to the cells at the edges, you should only count those on two sides to avoid the risk of overestimates or underestimates.

In the example shown above we perform the cell counting on a sample concentrated 10 times:

N = 5 (number of squares considered for counting)

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The device is manufactured under license Hycor Biomedical Inc.(EU Patent Number 326349)

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CLASSIFICATION : *In order to the Directive 98/79/CE concerning medical-diagnostic in vitro devices, Biosigma SpA, defines all the components of the above-state family as medical-diagnostic in vitro devices and classifies it as "Other medical-diagnostic in vitro devices", in order to make conformable the attachment III of the directive.*

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MATERIALS: *The materials utilized for the production of the in vitro labware devices are selected on the base of the properties indicated for the individual destination use.*

The device BVS100 is manufactured in METHACRYLATE.

Methacrylate is already known as raw material for manufacturing of in vitro medical device. Its mains characteristics are: rigid, transparent, resistant to atmospheric agents, it replace the glass in every its application in which it reach high temperature (lower than 90-100 °C).

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

The STAIN has the composition as follow:

| Denominazione | Conc. (C) | Denomination |
|--------------------------|----------------|-------------------------|
| VIOLETTO CRISTALLI | 0,05 ≤ C < 0,1 | VIOLET CRYSTALS |
| AMMONIO OSSALATO | 0 ≤ C < 0,05 | AMMONIUM OXALATE |
| SAFRANINA O | 0,24 ≤ C < 0,3 | SAFRANIN O |
| ETANOLO | 9,9 ≤ C < 11,4 | ETHANOL |
| INERTE(acqua distillata) | 87,9 ≤ C < 92 | INERT (distilled water) |

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The prepared is not classified as dangerous. The Security Data Sheet is available on request.

The reagent should be stored away from direct sunlight and at a temperature between +15°C and +25°C.

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IN VITRO DIAGNOSTIC MEDICAL DEVICE CE For professional use only

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RECOMMENDATIONS FOR USE (as for information)

- Use only for laboratory analysis.
- Do not hold the device to the flame or heat sources that may cause damage.
- Do not use after expiry or if box is already open.
- Do not vary the intended purpose of the product.
- Keep out of reach of children.
- Store in a dry place, away from sources of heat and direct sunlight.
- Waste disposal: utilize the provided D.P.I and discharge under the laws in force.

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The Declaration of Conformity in accordance to the Directive 98/79/CE is available, on request, to the Quality Department Biosigma SpA

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BIOSIGMA SpA - Via Valletta, 6 - 30010 - Cantarana di Cona (VE) - ITALY