Rev 24/01/2025

FLYFood

Bloomington Formulation: Powdered Drosophila Media by Biosigma

Description:

FLYFood Bloomington Formulation is meticulously designed to meet the needs of both large-scale Drosophila media preparation facilities and smaller individual laboratories. Packaged in a sturdy plastic bottle for easy handling and storage, this formulation guarantees consistent, high-quality food production. The recipe is versatile, allowing you to scale up (double, triple, or more) to meet any media quantity requirement. With a two-year shelf life from delivery, each batch is produced under strict quality control, ensuring rigorous ingredient traceability for unmatched reliability.

Ingredients: yeast, soy flour, yellow corn meal, agar, corn syrups solids*.

** Corn syrup solids replace light corn syrup.*

Instructions:

1. Pour the entire contents of the bottle (176g) into a suitable cooking vessel.
2. Gradually add ¼ or ½ of the total required amount of room temperature (RT) water (1 liter per bottle as reported in table 1) and stir to dissolve any clumps.
3. Add the remaining RT water and stir thoroughly.
4. Bring the mixture to a boil, then immediately reduce the heat and cover the vessel.
5. Let it simmer for 10 minutes, stirring frequently.
6. Once the mixture has cooled to approximately 70°C, add the specified amount of propionic acid (not supplied), as indicated in Table 1.
7. Stir well to ensure even distribution of the propionic acid.
8. Pour the prepared media immediately after adding the propionic acid.

Table1: Recommended Water and Propionic Acid Quantities per FLYFood Batch.

Amount of FLYFood	Total Amount Water	Amount of Propionic Acid**
1 bottle (176g)	1 liter	4.8 ml
2 bottles (352g)	2 liters	9.6 ml
5 bottles (880g)	5 liters	24 ml
10 bottles (1760g)	10 liters	48 ml

****** Use the specified amount of propionic acid only if using ≥99% propionic acid solution (13.4M).

Practical Suggestions:

- Dry food is often the source of issues. Small changes in the cooking process that lead to increased water loss can significantly impact food quality. Pay special attention to this, especially in dry climates.
- If you're new to FLYFood, start with small batches. Begin with one bottle before combining and using multiple bottles, to familiarize yourself with the process.
- To effectively remove clumps before and during cooking, use a whisk or hand mixer. A spoon can also work, but a whisk is most effective.
- The food is quite thick, even when hot. Ensure your pump can handle the consistency before preparing large batches.
- Tegosept may be used in addition to, but not as a substitute for, propionic acid. Propionic acid alone is considered sufficient as a preservative. Tegosept can slow Drosophila growth during the larval and pupal stages, so it should only be used when additional fungal protection is necessary or to address sporadic mold issues.
- To easily handle tegosept, prepare a solution by dissolving 25 g of the substance in 250 mL of ethanol. Add 10 mL of this solution per 1L of food.

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Biosigma S.p.A Dominique Dutscher Group