**Biosigma s.r.l.**

**a Dominique Dutscher Company**

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) |

Fax ++39 0426 302228 | SMS ++39 348 4077376 |

E-mail info@biosigmaeu.com | http://www.biosigma.com |

**M521D**

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10/02/2016

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**ENRICHMENT SYSTEM FOR THE PARASITOLOGY**

**„In Vitro“ Diagnostic Uses Only**

**KIT FOR 50 ANALISYS**

**Kit composition:**

REF	Description	Manufacturer	CE
<b>BR 102-B</b>	<b>ENRICHMENT SYSTEM:</b> - 50 sample vials with Bailenger solution - 50 filter vials	BioRepair GmbH Kirchenstr. 5, 74889 Sinsheim, Germany	<b>x</b>
<b>BSD858B</b>	<b>Medium B</b> Ethyle Acetate – 100 ml	<b>Titolchimica S.P.A</b> Pontecchio Pol. (RO) Italy	
<b>BSD858C</b>	<b>Medium C</b> Gram Lugol 1% - 15 ml	<b>Titolchimica S.P.A</b> Pontecchio Pol. (RO) Italy	<b>x</b>

**BAILENGER-METHOD**

**Instruction for use**

**Name and purpose of use of the enrichment system**

Optimized system for the concentration of parasitic elements in stool samples. The described system is based on the reliable MIF-Method (Sedimentation principle).

**Precautions**

The provided Bailenger solution contains acetic acid and sodium azide is toxic. In case of contact with skin immediately flush vigorously with water.

**Storage and expiry date**

Filled and empty vials are to be kept at room temperature until the expiry date on the packaging.

**Compounds of the enrichment system:**

1. filled sample vials with Bailenger-solution (Xi-irritant); (0,4 %; acetic acid; Sodium acetate 15 g/l, 0,05 % sodium azide; cap with spoon.
2. Filter vials for centrifugation with inserted filter

**Samples**

Samples should be fresh and untreated. Stool samples are supplied in normal commercial sample vials and taken from these for the parasitic enrichment.

**Sample preparation**

For optimal results take a spoon full of stool into a sample vial with 3.3 ml Bailenger-solution, See flow chart fig. 1 .

**Test performance**

Performance of the parasites enrichment from stool samples (parallel testing per sample increases the accuracy of the results)

- shake the closed vial for a short time ( (shaker, vortex or manually), to separate the sample from the spoon and to distribute the sample in the Bailenger-solution.
- Open the sample vial again and add 1.25 ml medium B (ethyl-acetate) (extraction of the residual lipid debris). Then close the vial again and shake the closed vial for a short time ( (shaker, vortex or manually),
- Remove the cap and screw the filter vial on to the sample vial and turn the coupled system upside down.
- centrifuge the sample at 1500 – 2000 x rpm for 5 – 10 minutes, unscrew the filter together with the sample vial and discharge them. Discharge supernatant. Proceed only with the centrifugation vial.
- take the sediment directly or with a drop of Lugol-staining-solution on a slide for microscopic examination.

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- *examine under the microscope for several times increases the accuracy of the examination, minimum two attempts !  
(Addition of saline could make the pipetting easier),*
- *Recommended Resolution: 100 x for screening and for larger objects: eggs and larvae.  
400-1000 x for Protozoa  
(fungi, and especially yeasts are a good indication for a sufficient enrichment of the sample)*

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Melvin, D. and M. M. Brooke. 1974. *Laboratory procedures for the diagnosis of intestinal parasites.* U.S. Department of Health, Education and Welfare, publication no. (CDC) 75-8282. Center for Disease Control, Atlanta, GA.

### Warnings and precautions:

*The enrichment system is designed for a mechanical enrichment of parasites in stool samples.  
Stool samples have always to be kept as infectious material. (Protecting wears , gloves).*

### Quality control

*Each lab should define its own quality control system. The application of an enrichment system is recommended in ring-trials.*

### Note for device combination:

*In the centrifuging step fitting accessories for the centrifuge have to be present. For the enrichment both fixed angles rotors and swinging rotors are suitable. In case of swinging rotors free swinging should be possible.*

### Notes for secure disposing:

*For disposing the enrichment system the regulations of the lab for disposing polluted waste should be observed.*

### Procedure for decontamination

*In case of unintended flow out of reagents the liquid should be wiped up and flushed vigorously with water.*

### General notes

- *This diagnostic tool and all its containing compound may only be used for scientific purposes or if notified only for in vitro diagnostic.*
- *The supplied Baillenger-solution contains acetic acid and sodium azide. Please observe appropriate safety measures. Avoid contact with skin!*
- *Never pipette with your mouth.*
- *Wear disposable gloves during the test procedure .*
- *Do not use or mix reagents from different lots.*
- *For the quality control the guidelines for medical laboratories have to be observed.*
- *The characteristic data of the kit as incubation times, incubation temperature and pipetting volumes have been defined internally by the manufacturer. Changes not confirmed by the manufacturer may influence the results. The manufacturer can therefore not be held reliable for any resulting damage from this.*
- *In case of claims to guaranty the rejected product has to be resend to the manufacturer together with a written explanation within 14 days.*
- *Do not the use the kit the after the expiry date printed on the packaging of the kit.*

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### Summary

#### Security

- sterile sample after loading, no cooling is required
- small amounts of organic solvents

#### Comfort


- Handling without objections
- closed system, Sample tube used as a cap
- no bad smell after loading
- the sample tubes

#### Quality

- standardization (defined filter pore size) and possible quality control
- optimised concentration of the parasitic elements !


Pic.1 Sampling

Pic.2 Loading the sample tube

Pic.3 Addition of Liquid B


Pic.4 Mixing

Pic.5 Filtration

Pic.6 Phase separation

#### BAILENGER

**! Attention irritant !** Do not breath vapours, wash contaminations vigorously with water.

#### Medium B:

Ethylacetat (Attention inflammable, keep away from open flames !)

<examin under the microscope