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Solution for Long-Term Freezing and Storage of Microbial Cultures

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Description: CLEARStable is an innovative system for the freezing and long-term preservation of bacterial and microbial cultures. The system is designed to eliminate the need for freeze-drying or frequent transfers, ensuring optimal preservation of the cultures. The system includes a cryobox containing 25 or 80 sterile vial and 1 vial picker. Each vial includes 20 ± 1 porous beads, immersed in 1 ± 0.1 ml of a glycerol-based cryoprotective solution. Each production lot is subjected to a sterility test, with a Certificate of Analysis (CoA) issued to verify compliance with quality standards.

Purpose: CLEARStable ensures optimal preservation of bacteria, bacterial spores, fungi, mold and yeasts at low temperatures, maintaining their viability and genetic stability. It is ideal for microbiology laboratories, research, and biotechnology production.

Advantages:

- Long-term preservation: Maintains cultures for extended periods without loss of viability.

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Materials:

- **Cryobox:** Container designed for low-temperature storage made of polycarbonate.
- **Cryotubes:** Tubes made from medical-grade polypropylene with bi-injected polyethylene/rubber caps.
- **Beads:** Porous ceramic.
- **Cryoprotectant:** Glycerol-based solution.

Features:

- **Storage:** Before use, CLEARStable can be stored at room temperature or, preferably, at +4°C for optimal product preservation.
- **Storage Conditions for Microbiological Cultures:** Microorganisms can be stored at different temperatures depending on their specific biological nature. For short periods, storage at -20°C is recommended. For long-term storage, the ideal temperature is -80°C, while for prolonged conservation, the gaseous phase of liquid nitrogen can be used. The selection of the appropriate temperature depends on the biological characteristics of the microorganisms and the operating procedures implemented by the laboratories to ensure maximum viability and integrity of the cultures.
- **Sterility Test:** Performed with incubation for 3 days at 30°C and 5 days at 25°C. Final check at 7 days. No microbial growth detected.
- **Performance test** performed on *Escherichia coli* XL1-Blue, *Escherichia coli* ATCC 8739, *Staphylococcus aureus* ATCC 6538, *Pseudomonas aeruginosa* ATCC 9027, *Bacillus subtilis* ATCC 6633, *Candida albicans* ATCC 10231, *Aspergillus brasiliensis* ATCC 16404, *Clostridium sporogenes* ATCC 11437. Incubation under temperature and duration conditions specific to each microorganism. Adequate microbial growth observed for all tested organisms.
- **Validity:** 24 months.

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