



Biosigma SpA

a Dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.)

| Fax ++39 0426 302228 | SMS ++39 348 4077376 |

E-mail info@biosigmaeu.com | <http://www.biosigma.com> |

M521D

Rev. 30/10/2023

SCHEDA TECNICA PRODOTTO / TECHNICAL DATA SHEET

REF: BSA021

Conetto tipo OLLY / Cones OLLY Type



PER USO GENERICO / FOR GENERIC USE

Caratteristiche		Features
Diametro	10 mm	<i>Diameter</i>
Altezza	40 mm	<i>Height</i>
Volume	1,5 ml	<i>Volume</i>
graduazioni	Si/Yes	<i>Graduations</i>
Materiale	PP	<i>Raw material</i>
Sterile	NO	<i>Sterile</i>
Prodotto IVD, marcato CE	Si/Yes	<i>IVD item, CE marked</i>
Conf. pezzi	10 X 1.000	<i>Conf. Pcs.</i>

DEFINIZIONE DEL FABBRICANTE/ MANUFACTURER DEFINITION

BIOSIGMA SpA

Via Valletta, 6 - 30010 - Cantarana di Cona (VE) - ITALY

Tel. +39.0426.302224 - Fax +39.0426.302228

E-mail: info@biosigma.com - <http://www.biosigma.com>

Destinato esclusivamente ad uso professionale

- In caso di fuoriuscita del contenuto usare i guanti: pericolo di contaminazione
- Non avvicinare il dispositivo alla fiamma o a fonti di calore che lo potrebbero danneggiare.
- Non utilizzare il prodotto scaduto
- Non variare la destinazione d'uso
- Prodotto non adatto ai bambini
- Conservare in luogo asciutto, temperatura min -10° max +50°
- Smaltimento: utilizzare gli appositi D.P.I e smaltire secondo la normativa vigente

For professional use only

- *Use gloves: contamination risk if contents leaks*
- *Keep out of flame or heat sources which might damage the product*
- *Do not use after expiry date or if packing is opened*
- *Do not vary the intended purpose of the product*
- *Keep out of reach of children*
- *Store in dry place, temperature range: min. - 10 °C max 50°C*
- *Disposal: use appropriate personal protective equipment and act according to applicable regulations*