



**Biosigma S.p.A.**

**a Dominique Dutscher Company**

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) |

Fax ++39 0426 302228 | SMS ++39 348 4077376 |

E-mail info@biosigmaeu.com | http://www.biosigma.com | www.biosigma.it

**M521D**

Rev. 06/07/2021

**SCHEDA TECNICA PRODOTTO / TECHNICAL DATA SHEET**

**SACCHETTI PER TRASPORTO CAMPIONI  
SPECIMEN TRANSPORT BAGS**



**SACCHETTI IN LDPE A DOPPIA TASCA PER IL TRASPORTO DEI CAMPIONI, TASCA SEPARATA PER I RELATIVI DOCUMENTI. CON SIMBOLO DI RISCHIO BIOLOGICO E AREA DI SCRITTURA.**

**CODICE CND W05020180**

**CONTENITORI PER TRASPORTO DI CAMPIONI - ACCESSORI**

**DOUBLE POCKET BAGS IN LDPE FOR THE TRANSPORT OF SAMPLES, SEPARATE POCKET FOR THE DOCUMENTS. WITH BIOHAZARD SYMBOL AND WRITING AREA.**

Prodotto con marchio CE - conforme alla Direttiva 98/79/CE e al D.lgs 332 del 08/09/2000  
CE Marked product - manufactured in compliance with 98/79/CE Directive



REF	Colore / Color	Conf.
BSA2401	Neutro / Neutral	1000
BSA2101	Giallo / yellow	1000
BSA2201	Verde / green	1000
BSA2301	Azzurro / light blue	1000
BSA2461	Rosso / red	1000
BSA2601	Arancio/orange	1000
BSA2701	Rosa/Pink	1000
BSA2801	Viola/Violet	1000

Per ridurre i notevoli rischi di infezione nella manipolazione e nel trasporto di campioni biologici è stato realizzato un nuovo tipo di sacchetto in polietilene a bassa densità costituito da due tasche appositamente studiate.

- La prima, completamente trasparente, con zip a tenuta, di formato cm 16 x16 spessore 50my (tolleranza +/- 5%), può contenere fino a 8 provette da 10 ml o fino 10 da 5 ml.
- La seconda, colorata, di formato cm 16x23 spessore 30my (tolleranza +/-5%), senza chiusura e saldata all'altra tasca, contiene il modulo di richiesta esami.

*To reduce the significant risks of infection in the handling and transport of biological samples has been developed a new type of low-density polyethylene bag consists of two special pockets.*

*The first, completely transparent, with zip, format 16 x 16 cm thickness 50my (tolerance +/- 5%), it can hold up to 8 tubes of 10 ml or up to 10 tubes of 5 ml.*

*The second, welded to the first, transparent or colored, of size 16x23 cm thickness 30my (tolerance +/- 5%), without closure to contain the required diagnostic analysis document.*

**DESCRIZIONE:**

- Busta in polietilene con cerniera richiudibile a pressione ottenuta in estrusione.
- Area di scrittura.
- Indicazioni sull'utilizzo del sacchetto e condizioni di conservazione del campione (congelamento, refrigerazione, a temperatura ambiente), stampate sul sacchetto.
- Presenza del simbolo di "Rischio Biologico" in nero, conforme al "Biological Risk" in accordo alla EN 980:2008 – Symbols for use in the labelling of medical devices.



**Biosigma S.p.A.**

**a Dominique Dutscher Company**

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) |

Fax ++39 0426 302228 | SMS ++39 348 4077376 |

E-mail [info@biosigmaeu.com](mailto:info@biosigmaeu.com) | <http://www.biosigma.com> | [www.biosigma.it](http://www.biosigma.it)


**M521D**

Rev. 06/07/2021

## SCHEDA TECNICA PRODOTTO / TECHNICAL DATA SHEET

### DESCRIPTION:

- Polyethylene bags with pressure resealable zipper obtained by extrusion.
- Writing area.
- Instructions on using the bag and sample storage condition legend (frozen, refrigerated, room temperature), printed on the bag.
- Presence of the symbol of "Biohazard", in black, in conformity with the "Biological Risk" according to EN 980: 2008 - Symbols for use in the labeling of medical devices.



### TOLLERANZE DIMENSIONALI

Le tolleranze dimensionali possono variare di +/- 5 %.

### DIMENSIONAL TOLERANCES

Dimensional tolerances may vary by +/- 5%.

### MATERIALE

I sacchetti sono realizzati LDPE (LATEX FREE) Il film è estruso.

Le schede di sicurezza del LDPE sono rilasciate dal fornitore e disponiamo di una copia delle stesse presso i nostri uffici.

### MATERIAL

Bags are made with a LDPE (LATEX FREE). The film is bubbled extruded.

The LDPE MSDS are issued by the supplier and we have a copy of them at our offices.

Prodotto con marchio CE - conforme alla Direttiva 98/79/CE e al D.lgs 332 del 08/09/2000

CE Marked product - manufactured in compliance with 98/79/CE Directive

**DEFINIZIONE DEL FABBRICANTE/  
MANUFACTURER DEFINITION:**

**BIO SIGMA SpA**  
Via Valletta, 6  
30010 - Cantarana di Cona (VE) - ITALY  
Tel. +39.0426.302224 - Fax +39.0426.302228  
E-mail: [info@biosigma.com](mailto:info@biosigma.com)  
<http://www.biosigma.com>