



Biosigma SpA

a Dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) |

Fax ++39 0426 302228 | SMS ++39 348 4077376 |

E-mail info@biosigmaeu.com | <http://www.biosigma.com> |

M521D

Ultimo aggiornamento
19/05/2023

SCHEDA TECNICA PRODOTTO / TECHNICAL DATA SHEET

BSP048

Provetta cilindrica 10 ml
Cylindrical test tube 10 ml
Conf. 2.000 pz/pcs.



BSP049

Provetta cilindrica 10 ml - **etichettata**
*Cylindrical test tube 10 ml - with **label***
Conf. 2.000 pz/pcs



Altezza provetta	100 mm	Tube height
Diametro provetta	16 mm	Tube diameter
Gradazione	NO	Scale area
tappata	NO	With cap
Materiale	PP	Raw material
Sterile	NO	Sterile
Confezionamento	2000 pz/pcs	Package

Destinato esclusivamente ad uso professionale

- **Per uso generico di laboratorio**
- In caso di fuoriuscita del contenuto usare i guanti: pericolo di contaminazione
- Non avvicinare il dispositivo alla fiamma o a fonti di calore che lo potrebbero danneggiare
- Non utilizzare il prodotto scaduto
- Non variare la destinazione d'uso
- Prodotto non adatto ai bambini
- Conservare in luogo asciutto, temperatura min - 10 °C max + 50 °C
- Smaltimento: utilizzare gli appositi D.P.I. e smaltire secondo la normativa vigente

For professional use only

- **For general laboratory use**
- Use gloves: contamination risk if contents leaks
- Keep out of flame or heat sources which might damage the product
- Do not use after expiry date
- Do not vary the intended purpose of the product
- Keep out of reach of children
- Store in dry place, temperature range: min. - 10°C max 50°C
- Disposal: use appropriate personal protective equipment and act according to applicable regulations

BIOSIGMA SpA

Via Valletta, 6

30010 - Cantarana di Cona (VE) - ITALY

Tel. +39.0426.302224 - Fax +39.0426.302228

E-mail: info@biosigma.com

www.biosigma.it ; www.biosigma.com