



Biosigma SpA

a Dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.)

| Fax ++39 0426 302228 | SMS ++39 348 4077376 |

E-mail info@biosigmaeu.com | <http://www.biosigma.com> |

M521D


Ultimo aggiornamento
31/10/2023

SCHEDA TECNICA PRODOTTO/ TECHNICAL DATA SHEET

REF: BSP054

Provetta cilindrica Vol. 10 ml, Ø 16 X 100 mm - in PMMA

Cylindrical Test tube Ø 16x100 mm, 10 ml in PMMA



Famiglia: PROVETTE e TAPPI

REF	Descrizione/Description	Dimens.	Conf.	Sterile	Vel. centrifug.
BSP054	Provetta cilindrica 10 ml in PMMA <i>Cylindrical Test tube 10 ml in PMMA</i>	Ø 16 x 100 mm	2000	no	2000 x g
BSP055	Provetta cilindrica 10 ml in PMMA con etichetta <i>Cylindrical Test tube 10 ml in PMMA with Label</i>	Ø 16 x 100 mm	2000	no	2000 x g

PER USO GENERICO

Destinato esclusivamente ad uso professionale

RACCOMANDAZIONI D'USO (a carattere esclusivamente informativo)

- In caso di fuoriuscita del contenuto usare i guanti pericolo di contaminazione
- Non utilizzare il prodotto scaduto
- Non variare la destinazione d'uso
- Prodotto non adatto ai bambini
- Conservare in luogo asciutto, lontano da fonti di calore e dalla luce solare diretta
- Smaltimento: utilizzare gli appositi D.P.I. e smaltire secondo la normativa vigente

FOR GENERIC USE

for professional use only

RECOMMENDATION USE (exclusive as information)

- Use gloves: contamination risk if contents leaks
- Do not use after expiry date
- Do not vary the intended purpose of the product
- Keep out of reach of children
- Store in dry place and far from sources of heat and from direct sunlight
- Waste disposal: use the provided D.P.I. and discharge under the laws in force.



BIOSIGMA SpA

Via Valletta, 6 - 30010 - Cantarana di Cona (VE) - ITALY

Tel. +39.0426.302224 - Fax +39.0426.302228

E-mail: info@biosigma.com - <http://www.biosigma.it> - <http://www.biosigma.com>