



Biosigma S.p.A.

a Dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) |

Fax ++39 0426 302228 | SMS ++39 348 4077376 |

E-mail info@biosigmaeu.com | <http://www.biosigma.com> |

M521D

Ultimo aggiornamento
01/06/2023

SCHEDA TECNICA PRODOTTO / TECHNICAL DATA SHEET

BSP060

Provetta conica 10 ml
Conical test tube 10 ml
Conf. 2.000 pz/pcs.

BSP061

Provetta conica 10 ml - etichettata
Conical test tube 10 ml - with label
Conf. 2.000 pz/pcs



Altezza provetta	103 mm	Tube height
Diametro provetta	16 mm	Tube diameter
Volume	10 ml	Volume
Gradazione	SI / YES	Scale area
tappata	NO	With cap
Materiale	PP	Raw material
Sterile	NO	Sterile
Confezionamento	2000 pz/pcs	Package

Destinato esclusivamente ad uso professionale

- **Per usi generici di laboratorio**
- In caso di fuoriuscita del contenuto usare i guanti: pericolo di contaminazione
- Non avvicinare il dispositivo alla fiamma o a fonti di calore che lo potrebbero danneggiare
- Non utilizzare il prodotto scaduto
- Non variare la destinazione d'uso
- Prodotto non adatto ai bambini
- Conservare in luogo asciutto, temperatura min - 10 °C max + 50 °C
- Smaltimento: utilizzare gli appositi D.P.I. e smaltire secondo la normativa vigente

for professional use only

- **For general labware use**
- Use gloves: contamination risk if contents leaks
- Keep out of flame or heat sources which might damage the product
- Do not use after expiry date
- Do not vary the intended purpose of the product
- Keep out of reach of children
- Store in dry place, temperature range: min. - 10°C max 50°C
- Disposal: use appropriate personal protective equipment and act according to applicable regulations

DEFINIZIONE DEL FABBRICANTE / MANUFACTURER DEFINITION

BIOSIGMA s.r.l.

Via Valletta, 6 – 30010 - Cantarana di Cona (VE) - ITALY

Tel. +39.0426.302224 - Fax +39.0426.302228

E-mail: info@biosigma.com - www.biosigma.it ; www.biosigma.com