



Biosigma SpA

a Dominique Dutscher Company

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) |

Fax ++39 0426 302228 | SMS ++39 348 4077376 |

E-mail info@biosigmaeu.com | http://www.biosigma.com |

M521D

Ultimo aggiornamento

Last updating

17/05/2022

SCHEMA TECNICA PRODOTTO - TECHNICAL DATA SHEET

REF: BSR001 (conf./pkg. 25 x 1.000 pz/pcs) – **BSR001R** (conf./pkg. 1.000 pz/pcs)

Puntali tipo Eppendorf 0 - 200 µl

Puntale giallo tipo Eppendorf-Brand-Socorex da 0 a 200 µl

Tips Eppendorf type 0 - 200 µl

Yellow pipet tip Eppendorf-Brand-Socorex type from 0 to 200 µl



Tipo Eppendorf	MODELLO/ MODEL	Eppendorf Type
DIMENSIONI ESTERNE		EXTERNAL DIMENSIONS
Capacità	200 µl	Capacity
Diametro inferiore	N.A.	Lower diameter
Diametro superiore	8.3 mm	Upper diameter
Altezza	49 mm	Height
Materiale	PP	Raw material
Colore	GIALLO / YELLOW	Color
Sterile	NO	Sterile
Prodotto IVD, marcato CE	NO	IVD item, CE marked
BSR001 Conf. pz.	25 X 1.000	Pkg. Pcs. BSR001
BSR001R Conf. pz.	1.000	Pkg. Pcs. BSR001R

RACCOMANDAZIONI D'USO (a carattere informativo)

- Per usi generici di laboratorio
- In caso di fuoriuscita del contenuto usare i guanti pericolo di contaminazione
- Non utilizzare il prodotto scaduto
- Conservare in luogo asciutto, lontano da fonti di calore a dalla luce solare diretta
- Smaltimento: utilizzare gli appositi D.P.I e smaltire secondo la normativa vigente

RECOMMENDATION USE (exclusive as information)

- For general uses
- Use gloves: contamination risk if contents leaks
- Do not use after expiry date
- Do not vary the intended purpose of the product
- Keep out of reach of children
- Store in dry place and far from the direct light
- Waste disposal: use the provided D.P.I and discharge under the laws in force.