



**Biosigma SpA**

**a Dominique Dutscher Company**

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) |  
Fax ++39 0426 302228 | SMS ++39 348 4077376 |  
E-mail info@biosigmaeu.com | http://www.biosigma.com |

**M521D**

Ultimo aggiornamento  
Last updating  
17/05/2022

## SCHEMA TECNICA PRODOTTO - TECHNICAL DATA SHEET

**REF: BSR002** (conf./pkg. 20 x 500 pz/pcs) – **BSR002R** (conf./pkg. 500 pz/pcs)

**Puntali tipo Eppendorf 100 - 1000 µl**

Puntale azzurro tipo Eppendorf-Brand-Socorex  
da 200 a 1000 µl

**Tips Eppendorf type 100 - 1000 µl**

Light blue pipet tip Eppendorf-Brand-Socorex type  
from 200 to 1000 µl



<b>Tipo Eppendorf</b>	<b>MODELLO/ MODEL</b>	<b>Eppendorf Type</b>
<b>DIMENSIONI ESTERNE</b>		<b>EXTERNAL DIMENSIONS</b>
Capacità	<b>1000 µl</b>	Capacity
Diametro inferiore	<b>N.A.</b>	Lower diameter
Diametro superiore	<b>9.5 mm</b>	Upper diameter
Altezza	<b>71 mm</b>	Height
Materiale	<b>PP</b>	Raw material
Colore	<b>Azzurro / Light blue</b>	Color
Sterile	<b>NO</b>	Sterile
Prodotto IVD, marcato CE	<b>NO</b>	IVD item, CE marked
BSR002 Conf. pz.	<b>20 X 500</b>	Pkg. Pcs. BSR002
BSR002R Conf. pz.	<b>500</b>	Pkg. Pcs. BSR002R

**RACCOMANDAZIONI D'USO (a carattere informativo)**

- Per usi generici di laboratorio
- In caso di fuoriuscita del contenuto usare i guanti pericolo di contaminazione
- Non utilizzare il prodotto scaduto
- Conservare in luogo asciutto, lontano da fonti di calore a dalla luce solare diretta
- Smaltimento: utilizzare gli appositi D.P.I e smaltire secondo la normativa vigente

**RECOMMENDATION USE (exclusive as information)**

- For general uses
- Use gloves: contamination risk if contents leak
- Do not use after expiry date
- Do not vary the intended purpose of the product
- Keep out of reach of children
- Store in dry place and far from the direct light
- Waste disposal: use the provided D.P.I and discharge under the laws in force.