



**Biosigma SpA**

**a Dominique Dutscher Company**

Via Valletta, 6 | 30010 Cantarana di Cona (VE), Italy | Tel. ++39 0426 302224 (r.a.) |

Fax ++39 0426 302228 | SMS ++39 348 4077376 |

E-mail info@biosigmaeu.com | http://www.biosigma.com |

**M521D**

Ultimo aggiornamento

Last updating

17/05/2022

## SCHEMA TECNICA PRODOTTO - TECHNICAL DATA SHEET

**REF: BSR066** (Q.tà/Q.ty 25 x 1000 pz./pcs.) – **BSR066R** (Q.tà/Q.ty 1000 pz./pcs.)



**Puntali gialli tipo Gilson 0 - 200 µl**

**Yellow tips Gilson type 0 - 200 µl**

Tipo Gilson	MODELLO/ MODEL	Gilson Type
DIMENSIONI ESTERNE		EXTERNAL DIMENSIONS
Capacità	200 µl	Capacity
Diametro inferiore	N.A.	Lower diameter
Diametro superiore	7 mm	Upper diameter
Altezza	50 mm	Height
Materiale	PP	Raw material
Colore	Giallo / Yellow	Color
Sterile	NO	Sterile
Prodotto IVD, marcato CE	NO	IVD item, CE marked
BSR066 Quantità pz.	25 X 1.000	Quantity pcs. BSR066
BSR066R Quantità pz.	1.000	Quantity pcs. BSR066R

### RACCOMANDAZIONI D'USO (a carattere informativo)

- Per usi generici di laboratorio
- In caso di fuoriuscita del contenuto usare i guanti pericolo di contaminazione
- Non utilizzare il prodotto scaduto
- Conservare in luogo asciutto, lontano da fonti di calore e dalla luce solare diretta
- Smaltimento: utilizzare gli appositi D.P.I e smaltire secondo la normativa vigente

### RECOMMENDATION USE (exclusive as information)

- For general uses
- Use gloves: contamination risk if contents leak
- Do not use after expiry date
- Do not vary the intended purpose of the product
- Keep out of reach of children
- Store in dry place and far from the direct light
- Waste disposal: use the provided D.P.I and discharge under the laws in force.