





BIODipper SAMPLE CONTAINERS

Single-use immersion containers for sampling and secure transport of the sample in the same unit, eliminating the risk of cross-contamination.

BIODipper Sample Containers are single-use immersion containers designed to streamline the sampling process, enhancing both safety and efficiency.

They allow the collection, safe transport, and direct analysis of the sample using a single container, significantly reducing the risk of cross-contamination during laboratory procedures.


A practical and reliable solution, ideal for ensuring sample integrity throughout every stage of the analysis.





Characteristics:

- **Transparent blue color**, perfect for use in the food, feed and pharmaceutical industries.
- Ideal for **Beverage Manufacturing Industry**
- Suitable for food contact.
- The **ergonomic handle** is designed with a **grip ring** for safe and secure handling.
- 5 available volumes: 20mL, 50mL, 60mL, 120mL e 160mL.
- 2 handle lengths available: 25cm e 35cm.
- Pack: single.
- After sampling, **the container is sealed**, the **locking strap and handle can be removed** to simplify transport and analysis, and the sample can be safely sent to the laboratory.
- Treated by **Gamma rays** irradiation at 25kGy, sterility test performed on each production batch.
- Validity: 48 months.
- Ideal for extracting samples from barrels, tanks and other containers.



ITEM	DESCRIPTION	HANDLE	GRADUATION	PACK
DIP0020L	Sampling container 20 ml, handle in PA, blue	short	20mL	50
DIP0020H	Sampling container 20 ml, handle in PA, blue	long	20mL	50
DIP0050L	Sampling container 50 ml, handle in PA, blue	short	50mL	50
DIP0050H	Sampling container 50 ml, handle in PA, blue	long	50mL	50
DIP0060L	Sampling container 60 ml, handle in PA, blue	short	60mL	40
DIP0060H	Sampling container 60 ml, handle in PA, blue	long	60mL	40
DIP0120L	Sampling container 120 ml, handle in PA, blue	short	120mL	80
DIP0120H	Sampling container 120 ml, handle in PA, blue	long	120mL	80
DIP0160L	Sampling container 160 ml, handle in PA, blue	short	160mL	80
DIP0160H	Sampling container 160 ml, handle in PA, blue	long	160mL	80