



Biosigma

Deutscher group




STERILE INOCULATION LOOPS

Three loop types. Precision, control, and the right flexibility.

Made of high-quality materials and available in two volumes (1 μ L and 10 μ L), our inoculation loops feature a double-ended design – loop head and needle – ensuring maximum operational flexibility.

The easy-to-open peel-pack allows quick access even when wearing gloves, while lot number and expiry date guarantee full traceability and safe use in the lab.



STANDARD

Material: Polypropylene
Colors: white, blue

Balanced rigidity and ease of handling for routine microbiology. Ensures reliable, consistent streaking in daily lab procedures.

SOFT

Material: Polypropylene
Colors: light green, light blue

Flexible design for delicate media or tasks requiring gentle handling. Minimizes agar damage and provides excellent tactile feedback.

RIGID

Material: Polystyrene
Colors: green, violet

Firm structure for stable pressure and precise control. Ideal for demanding applications where consistency is critical.

Standard loops



Item	Color	Volume	Packaging	Sales unit
BSV121	White	1µL	peel pack, 20pcs	10.000
BSV1210	White	1µL	peel pack, 10pcs	12.000
BSV1215	White	1µL	peel pack, 5pcs	4.000
BSV1211	White	1µL	individually wrapped	2.100
BSV120	Blue	10µL	peel pack, 20pcs	10.000
BSV1200	Blue	10µL	peel pack, 10pcs	12.000
BSV1205	Blue	10µL	peel pack, 5pcs	4.000
BSV1201	Blue	10µL	individually wrapped	2.100

Soft loops



Item	Color	Volume	Packaging	Sales unit
BSV121S	Light green	1µl	peel pack, 20pcs	10.000
BSV1210S	Light green	1µl	peel pack, 10pcs	12.000
BSV1215S	Light green	1µl	peel pack, 5pcs	4.000
BSV1211S	Light green	1µl	individually wrapped	2.100
BSV120S	Light blue	10µl	peel pack, 20pcs	10.000
BSV1200S	Light blue	10µl	peel pack, 10pcs	12.000
BSV1205S	Light blue	10µl	peel pack, 5pcs	4.000
BSV1201S	Light blue	10µl	individually wrapped	2.100

Rigid loops



Item	Color	Volume	Packaging	Sales unit
BSV121R	Green	1µl	peel pack, 20pcs	10.000
BSV1210R	Green	1µl	peel pack, 10pcs	12.000
BSV1215R	Green	1µl	peel pack, 5pcs	4.000
BSV1211R	Green	1µl	individually wrapped	2.100
BSV120R	Violet	10µl	peel pack, 20pcs	10.000
BSV1200R	Violet	10µl	peel pack, 10pcs	12.000
BSV1205R	Violet	10µl	peel pack, 5pcs	4.000
BSV1201R	Violet	10µl	individually wrapped	2.100



Biosigma
Dutscher group

Biosigma S.p.A Dominique Dutscher Group
Via Valletta, 6 | 30010, Cona (Venice) Italy
Phone +39 0426 302226 | info@biosigmaeu.com
www.biosigma.it | www.biosigma.com

Rev 10/2025