
False Bottom Tubes

2.5 ml

Biosigma False Bottom Tubes are designed for low-volume sample handling, combining efficient recovery, reliable identification and flexible configurations to support your laboratory workflows.



CLEAR VERSION

Clarity for everyday laboratory workflows.

Designed for efficient handling of low-volume samples in routine laboratory applications.

AMBER VERSION

For light-sensitive samples.

An amber colored solution for laboratory applications involving light-sensitive samples.

CONICAL FALSE BOTTOM

Designed to concentrate small sample volumes in a reduced collection area.

2.5 ml	100 µl	13x75 mm
NOMINAL VOLUME	MINIMUM DEAD VOLUME	TUBE FORMAT

Recover more for every sample

The conical false-bottom design concentrates small sample volumes in a reduced collection area, supporting efficient sample handling and recovery. Ideal for low-volume and valuable samples.

Made for traceable workflows

A generous labelling area allows the application of labels and barcodes, supporting clear sample identification throughout laboratory workflows, from sample receipt and analysis to storage and archiving.

Designed around your laboratory

- Made from polypropylene (PP)
- 13 × 75 mm tube format with graduation
- Compatible with laboratory analysers and pre- and post-analytical automation
- Rounded outer base for stable positioning in racks
- Available in clear and amber versions
- Screw cap supplied separately, or pre-assembled

1 CHOOSE YOUR TUBE

Two versions for different laboratory needs.



CLEAR VERSION

For routine laboratory applications.



AMBER VERSION

For light-sensitive samples.



2 CHOOSE YOUR CLOSURE

Compatible screw caps, supplied separately or pre-assembled.



Neutral cap



Amber cap



Screw cap supplied separately or pre-assembled.

3 IDENTIFY THE SAMPLE

Colored cap disks for easy sample identification.



Orange | Light-blue | White | Yellow | Pink | Red | Green | Purple

ORDERING INFORMATION

False Bottom Tubes Ø 13 mm and compatible accessories

False Bottom Tubes - 2.5 ml

Code	Colour	Cap	Units/Box	Sales Units
BSN250	Neutral	Assembled	500	1000
BSN250AM	Amber	Assembled	500	1000
BSN2500	Neutral	Without cap	1000	2000
BSN2500AM	Amber	Without cap	1000	2000

Compatible Screw Caps - Ø 13 mm

Code	Colour	Units/Box	Sales Units
BSN13	Neutral	1000	5000
BSN13AM	Amber	1000	5000

Colored Cap Disks

Code	Colour	Sales Units
BSN001AR	Orange	1000 pcs
BSN001AZ	Light Blue	1000 pcs
BSN001BI	White	1000 pcs
BSN001G	Yellow	1000 pcs
BSN001RA	Pink	1000 pcs
BSN001RO	Red	1000 pcs
BSN001V	Green	1000 pcs
BSN001VI	Purple	1000 pcs